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PK7669303 10/2023

## Instructions for Use

**Items needed for application (not included):**

Trimmers or Clippers  
Soap and water wipes  
Drying towel

**START HERE:** Wash your hands before and after this procedure. If placing or removing this product on another person, wear gloves. Please read instructions, warnings, precautions, and recommendations entirely before placing the device.

**Indications for use:**

The device is intended for non-invasive urine output management, such as urinary incontinence, in users with male anatomy.

**Contraindications:**

- Users with urinary retention.

**Warnings:**

- Do not use on irritated skin. Examples include, but are not limited to, rashes, skin tears, or blisters.
- Do not cover fresh surgical wounds with the device.
- Do not use the device with any material or in any position that blocks the vents on the front of the device or does not allow for airflow through the device. Examples include, but are not limited to, tight clothing or briefs.
- To avoid potential skin injury, never pull the device directly away from the user. Always peel in the direction from head to foot.
- Stop use if an allergic reaction occurs.
- After use, this product may be a potential biohazard. Dispose of in accordance with applicable local, state, and federal laws and regulations.


**Precautions:**


- Not recommended for users who:
  - Are agitated, combative, or uncooperative and might remove the device.
  - Have frequent episodes of bowel incontinence without a fecal management system in place.
  - Have skin irritation or breakdown at the site.
- Do not use barrier creams, lotions, or ointments on the penis or pubic skin around the base of the penis when using the device. These may impede suction or weaken adhesive.
- Proceed with caution in users who have had recent surgery of the external male anatomy.
- Always check skin for irritation or breakdown and clean and dry the skin prior to placement of a new device.

**Recommendations:**

- Wash hands thoroughly before device placement.
- Ensure the device remains connected after turning the user, monitoring for pulling and tension on the device.
- Remove the device prior to walking.
- Check device placement and user's skin at least every 2 hours.
- Refer to the PureWick™ Urine Collection System IFU for instructions on suction tubing replacement.

SYMBOL	SYMBOL TITLE (STANDARD REFERENCE)	EXPLANATORY TEXT
Symbols immediately below are derived from the following standard: ISO 15223-1:2021 - Medical devices — Symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements.		
	Manufacturer (5.1.1)	Indicates the medical device manufacturer
	Lot number (5.1.5)	Indicates the manufacturer's batch code so that the batch or lot can be identified
	Catalog number (5.1.6)	Indicates the manufacturer's catalog number so that the medical device can be identified
	Use by (5.1.4)	Indicates the date after which the medical device is not to be used
	Non-sterile (5.2.7)	Indicates a medical device that has not been subjected to a sterilization process
	Do not use if package is damaged and consult instructions for use (5.2.8)	Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information
	Single use (5.4.2.)	Indicates a medical device that is intended for one single use only
	Consult instructions for use or consult electronic instructions for use (5.4.3)	Indicates the need for the user to consult the instructions for use
Symbols immediately below are derived from the following standard: ISO 7000 – Graphical symbols for use on equipment		
	Packaging unit (2794)	To indicate the number of pieces in the package
<b>Other</b>		
	Not made with Natural Rubber Latex	Indicates that the medical device is not made with natural rubber latex

**PureWick™ Urine Collection System Setup**

1. Please consult the PureWick™ Urine Collection System Instructions for Use as needed. Confirm the PureWick™ Urine Collection System, collection canister, and tubing are set up correctly. Turn the unit ON, using the ON/OFF switch.  
**NOTE:** Ensure all connections are air-tight and check for possible cracks, leaks, kinks, or blockages.
2. Securely connect the drainage fitting of the PureWick™ Male External Catheter to the open end of the collector tubing on the PureWick™ Urine Collection System.

**Placement of PureWick™ Male External Catheter**

3. Trim or clip the pubic hair to ensure the product securely adheres to the skin. Check the skin and do not use if there is irritation or breakdown. Clean the penis and the skin around the base of the penis with soap and water or soap and water wipes. Dry skin prior to positioning the device.  
**NOTE:** Moisture or soap residue on skin will weaken the adhesive.
4. Position the device, aligning drainage fitting towards the feet of the user. Slide the opening of the device over the penis until close to, but not touching, the skin around the base of the penis. Center penis within the opening. Do not place scrotum inside the device. Remove the bottom adhesive peel-off and press the device against the skin below the base of the penis. Remove the top adhesive peel-off and press the device against the skin above the base of the penis. Ensure device has fully adhered around the base of the penis by pressing down on the adhesive.

**Removal of PureWick™ Male External Catheter**

5. Gently lift the top edge of the adhesive base off the skin. In a slow, downward peeling motion, remove the device. If necessary, use a warm, wet cloth or pad to help loosen adhesive.  
**Warning:** To avoid potential skin injury, never pull the device directly away from the user. Always peel in the direction from head to foot.  
After use, this product may be a potential biohazard. Dispose of in accordance with applicable local, state, and federal laws and regulations.

**When to Replace PureWick™ Male External Catheter**

6. Replace the device at least every 24 hours or if soiled with feces, blood, or semen.

**Wall Suction**

This product is compatible with wall suction in hospital-like settings. If using this set-up, connect the canister to wall suction and set to a minimum of 40 mmHg continuous suction. Always use the minimum amount of suction necessary.

BD, the BD Logo and PureWick are trademarks of Becton, Dickinson and Company, or its affiliate. © 2023 BD and its subsidiaries. All rights reserved.

Manufacturer:  
C. R. Bard, Inc.  
8195 Industrial Boulevard  
Covington, GA 30014 USA  
1.844.823.5433  
bd.com

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